












Accessories for packaging

Pack type	Material	Closure	Simple accessories	Advanced accessories	Accessory info
5 L Jerrycan	PE	51 mm	Tap DIN 51 Cat.N° 29548.042		
5 L Jerrycan	PE	42 mm	TAP Cat.N° 29048.098 (by 10 units)		
10 L Jerrycan	PE	44 mm	Cap with tap for 10 L PE container (Technisolv type) Cat.N° 223-0059		
10 L Jerrycan	PE	51 mm	Tap DIN 51 Cat.N° 29548.042	OTAL hand Cat.N° 223-1112 or foot-pump Cat.N° 223-1122	
20 L Jerrycan	PE	61 mm	Tap DIN 61 Cat.N° 29548.043	OTAL hand or foot-pump	
25 L Jerrycan	PE	61 mm	Tap DIN 61 Cat.N° 29548.043	OTAL hand or foot-pump	
25 L drum	Steel	2" & 3/4"	Tap 29548.096 + Tri-Sure (fine thread) adapter Cat.N° 223-1622	SS Pump Cat.N° 223-1313 + barrel joint Cat.N° 223-1249 + antistatic kit Cat.N° 223-1391	
25 L drum with liner	Steel	2" & 3/4"	Tap 29548.096 + Mauser adapter Cat.N° 223-1621	PP pump Cat.N° 223-1220 + barrel joint Cat.N° 223-1262	
60 L drum	PE	2" & 3/4"	Tap 29548.096 + Mauser adapter Cat.N° 223-1621	PP pump Cat.N° 223-1221 + barrel joint Cat.N° 223-1262	
			Tap 29548.096 + Tri-Sure (coarse thread) adapter Cat.N° 223-1620	PP pump Cat.N° 223-1221 + barrel joint Cat.N° 223-1263	
60 L drum with liner	Steel	2"	Tap 29548.096 + Mauser adapter Cat.N° 223-1621	PP pump Cat.N° 223-1221 + barrel joint Cat.N° 223-1262	
200 L drum with liner	Steel	2"	Tap 29548.096 + Mauser adapter Cat.N° 223-1621	PP pump Cat.N° 223-1223 + barrel joint Cat.N° 223-1262	
200 L drum	PE	2"	Tap 29548.096 + Mauser adapter Cat.N° 223-1621	PP pump Cat.N° 223-1223 + barrel joint Cat.N° 223-1262	
200 L drum	Steel	2" & 3/4"	Tap 29548.096 + Tri-Sure (fine thread) adapter Cat.N° 223-1622	SS pump Cat.N° 223-1235 + barrel joint Cat.N° 223-1249 + antistatic kit Cat.N° 223-1391	

OTAL is a registered trademark from Bürkle GmbH, Tri-Sure is a registered trademark from CDS



Other accessories

Type	Closure	Accessories	Accessories Info
Taps	3/4"	HDPE white with thread Cat.N° 29548.096	
	2"	HDPE white with thread Cat.N° 29548.097	
	3/4"	Tap brass modified design Cat.N° 331603M	
	3/4"	Tap polypropylene Cat.N° 331582B	
	51 mm & 61 mm	Cap with pouring flex : (DIN 51) Cat.N° 29548.075 (DIN 61) Cat.N° 29548.086	
Drum keys for screwing/unscrewing			
Plastic drum keys	51 mm	For 10 L (DIN 51) plastic drums, 1 piece Cat.N° 29548.064	
	61 mm	For 20 and 25 L (DIN 61) plastic drums, 1 piece Cat.N° 29548.053	
Metal drum key	2" & 3/4"	For seals and 2" and 3/4" closures on metal drums, 1 piece Cat.N° 29548.061	
Glass rods		For opening ConvoL NORMADOSE ampoules, 10 pieces/pack Cat.N° 238-0000	
Adapter	45 mm	Adapter for dispenser S40 to GL45, PP Cat. N° 217-9271	
Tape		Tape PTFE Cat. N° 332362S	
Bottle carrier		For 2.5 L bottle Cat.N° 215-2734	
Pumps (piston)		Pump for solvents For 20/25 L containers Cat.N° 181-0303 For 205 L containers Cat.N° 181-0304	
		Pump for chemicals For 20/25 L containers Cat.N° 181-0301 For 205 L containers Cat.N° 181-0302	
		Adapter for 205 L pumps 2" Tri-Sure Cat.N° 181-0300 2" Mauser Cat.N° 224-0000	